## TECHNICAL DATA

	SMART	SMART <sup>s</sup>
Seat total width in mm	320 / 340 / 360 / 380 / 400 / 420 / 440 / 460 / 480 / 500 / 520	320 / 340 / 360 / 380 / 400 / 420 / 440 / 460 / 480 / 500 / 520
Road-ready width in mm	500 / 520 / 540 / 560 / 580 / 600 / 620 / 640 / 660 / 680 / 700	500 / 520 / 540 / 560 / 580 / 600 / 620 / 640 / 660 / 680 / 700
Folded width in mm	280 at 0° camber or 330 at 3° camber	280 at 0° camber or 330 at 3° camber
Seat depth in mm	380 / 400 / 420 / 440 / 460 / 480 / 500	380 / 400 / 420 / 440 / 460 / 480 / 500
Front seat height in mm	440 – 530	400 – 540
Rear seat height in mm	400 – 530	370 – 530
Overall height in mm	Rear seat height + back height	Rear seat height + back height
Backrest height in mm	Variable from 250 to 500, + 25	Variable from 250 to 500, + 25
Armrest height in mm	approx. 170 – 245	approx. 170 – 245
Length with / without footrest in mm	880 - 1,030 (depending on frame length, wheel size and axis position)	880 - 1,030 (depending on frame length, wheel size and axis position)
Userweight in kg	150	150
Chair weight in kg	From 9.5	From 10.5
Allowed total weight in kg	From 159.5	From 160.5
Transport weight in kg	7.4	6.5
Drive wheel in mm (inches)	600 (24") 635 (25") 660 (26")	560 (22") 600 (24") 635 (25") 660 (26")

Dimensions in mm, weights in kg, wheels in inches, speed in km/h, dimension tolerances  $\pm$  10 mm, Data subject to design changes.

## FRAME COLOURS



#### **MEYRA GmbH**

- Meyra-Ring 2
  D-32689 Kalletal-Kalldorf
- Tel +49 5733 922 0 Fax +49 5733 922 - 9311
- 🖂 info@meyra.de

www.meyra.de





# **PERFECT HANDLING** LIGHT WEIGHT



Unique side guards with height adjustable armrest for one-hand operation



The new pressure brake offers extremely low actuation forces



Perfect chassis setting thanks to individually adjustable Castor housing

3



FrogLegs in 5 colours, anodised aluminium alloy

### SMART<sup>F</sup> AND SMART<sup>S</sup> – THE NEW ACTIVE SEGMENT

The new foldable adaptive wheelchairs with integrated leg support are the highlight of the active section and combine lightness with perfect handling to be the perfect companions in all life situations. Perfectly adaptable to physical changes and the individual needs of the user, comibined with perfect suspension settings. Design and functionality combined in an impressive wheelchair!





Magnesium wheels in 5 colours, painted hub and rim

Optimised folding system with low folding dimension offers a nearly parallel folding kinematic, enables independent installation

Function elements with one-hand operation newly designed side part with heightadjustable armrest

Individuality through colour concepts new design selection for Castor and drive wheels

Newly designed rotating leg supports and aluminium foot plates angle and depth-adjustable leg support, light and extremely stable in high quality

