

NOT ONLY FOR HOME...

... also the ideal companion on the way!

Thanks to compact dimensions and easy and fast dismantling is the MICRO-ELECTRO-MOBIL the faithful companion for on the go. Enjoy a trip to the shopping mall – without any problems when walking.

PERFECT FOR HOLIDAYS IN THE HOTEL



TECHNICAL DATA

	MICRO-ELECTRO-MOBIL / MICRO-ELECTRO-MOBIL UP
Length in mm	760
Width in mm	580
Adjustable seat height in mm	370 – 450
Seat depth in mm	400
Back height in mm	400
Adjustable seat width in mm	470 – 570
Transport height (with folded steering column and without seat) in mm	470
Number of wheels	5
Steering wheel diameter in mm	110
Drive wheel diameter in mm	220
Ground clearance in mm	50
Obstacle avoidance max. in mm	15
Turning circle (radius) in mm	500
Permissible gradients / slope in degrees (percent)	3° (6%)

	MICRO-ELECTRO-MOBIL	MICRO-ELECTRO-MOBIL UP
Continuous motor power in W	200	
Peak motor power in W	800	
Batteries in Ah	20,8Ah with lithium batteries 14Ah with lead batteries	
Battery charger in A	2	
Electronics in A	50	
Range, max. (under test conditions) in km	17 with lithium batteries 10 with lead batteries	
Tare weight in kg with battery	42 with lithium batteries 47 with lead batteries	47 with lithium batteries 52 with lead batteries
Maximum user weight in kg	120	
Maximum total weight in kg	170	
Speed, max. in km/h	6	
Lift height in mm	-	100

MEYRA GmbH

Meyra-Ring 2
D-32689 Kalletal-Kalldorf

Tel +49 5733 922 - 0
Fax +49 5733 922 - 9311

info@meyra.de

www.meyra.de

AMG ACTIVE MOBILITY GENERATION

The MEYRA MICRO-ELECTRO-MOBIL is part of the AMG concept.



MICRO-ELECTRO-MOBIL
The space saver



iTravel
The ideal travelling companion



iMove
The active footrest

The MEYRA series for independent living within your own four walls.



by www.bwcmgmbh



A REAL SPACE SAVER

The agile MICRO-ELECTRO-MOBIL for smallest rooms

The MICRO-ELECTRO-MOBIL is designed for people with reduced mobility. Thanks to compact dimensions and by means of Mid wheel drive, the MICRO-ELECTRO-MOBIL has a turning radius of 500 mm and is therefore perfectly suitable for indoor use.



FACILITATES EVERYDAY LIFE

Manoeuvrable

It turns on the spot by the small turning radius of 500 mm. Thus, small rooms and narrow doors are no problem.

Mid wheel drive

The Mid wheel drive is a system of class. Also carpets with a deep pile are no obstacle.

Turnable

With the turnable and height-adjustable seat, domestic work can be done easily in a sitting position.

Folding & compact

Through the simple and fast dismantling in three parts, it can be easily loaded, also in small cars.

Strong

Powerful 200 Watt engines allow a user weight of 120 kg.

Easy and simple

The operation is self-explanatory and so easy to understand. Easy recharging overnight guaranteed driving pleasure every day.

WITH LIFT: MICRO-ELECTRO-MOBIL UP

The lift offers the possibility to carry out various activities in the living environment, e.g. cooking, morning hygiene, reaching higher cabinets or setting the perfect sitting position at tables.